



SLS Lab Pro Pipettes

Single Channel

- Ergonomic design ensuring light weight & soft plunger movement
- Click-stop system for volume selection
- Colour coded for easy identification
- Calibration confirms to EN-ISO 8655 with report enclosed with every pipette

| Code | Channel | Volume Range, μL | Increment, μL |
|---------|---------|-----------------------------|--------------------------|
| SLS4900 | Single | 0.1-2.5 | 0.01 |
| SLS4902 | Single | 0.5-10 | 0.1 |
| SLS4910 | Single | 2-20 | 0.1 |
| SLS4904 | Single | 5-50 | 0.5 |
| SLS4906 | Single | 10-100 | 0.5 |
| SLS4912 | Single | 20-200 | 1 |
| SLS4908 | Single | 100-1000 | 5 |
| SLS4914 | Single | 500-5000 | 50 |
| SLS4916 | Single | 1000-10000 | 100 |



Multichannel

- Smooth and light tip loading and ejection
- Leak free sealing of tips
- Manifold can rotate 360° for easy right or left hand operation
- Easily removable & fully autoclavable manifold

| Code | Channel | Volume Range, μL |
|---------|---------|-----------------------------|
| SLS4926 | 8 | 0.5 - 10 |
| SLS4918 | 8 | 5 - 50 |
| SLS4920 | 8 | 10 - 100 |
| SLS4922 | 8 | 20 - 200 |
| SLS4924 | 8 | 30 - 300 |
| SLS4936 | 12 | 0.5 - 10 |
| SLS4928 | 12 | 5 - 50 |
| SLS4930 | 12 | 10 - 100 |
| SLS4932 | 12 | 20 - 200 |
| SLS4934 | 12 | 30 - 300 |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.