

1. Part Number: A643-500

2. Photograph:



3. Product Information:

Description:	Industrial Scoop
Nominal Scoop Volume:	500ml
Material of Construction:	Stainless Steel
Method of Construction:	Fully welded handle
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	290mm
Body Length:	185mm
Max. Body Width:	105mm
Max. Body Height:	65mm
Nominal Weight:	235g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jan-20



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.