


1. Part Number: 1070A-350

2. Photograph:



3. Product Information:

Description	MicroMate 350mm
Materials of Construction:	316L stainless steel Nylon Clutch (non product contact part)
Nominal Sample Volume:	Variable up to 0.1ml
Overall Length:	500mm
Length of Sampling Tube:	350mm
Diameter of Sampling Tube:	5mm
Weight:	200g
Finish:	<0.8 microns Ra
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.