

GyroStir Magnetic Stirrers & Hotplate Stirrers



SciQuip GyroStir Basic – Analogue 5” Magnetic Stirrer

- ABS acid resistant plastic 5” plate
- Speed range 0 – 1,500 rpm
- Max stirring capacity 3 litres (water)
- Classic analogue scale control
- Supplied with 1x stirring bar



Ordering Information

Part No.	Description
SP2201-B-A	GyroStir Basic - Analogue 5” Magnetic Stirrer

Technical Specifications

Workplate Dimensions	Ø 135 mm
Max Stirring volume (H2O)	3L
Speed Range	0 – 1,500 rpm
Display	Analogue Scale
Connectivity port	No
Dimensions	150 x 260 x 80 mm
Power	200-240V / 50-60Hz / 10W

SciQuip GyroStir Digital Basic – Digital 5" Magnetic Stirrer

- ABS acid resistant plastic 5" plate
- Speed range 100 – 1,500 rpm
- Max stirring capacity 3 litres (water)
- LED display – both set & real time values
- Supplied with 1x stirring bar



Ordering Information

Part No.	Description
SP2201-B-D	GyroStir Digital Basic - 5" Magnetic Stirrer

Technical Specifications

Workplate Dimensions	Ø 135 mm
Max Stirring volume (H2O)	3L
Speed Range	100 – 1,500 rpm
Display	LED
Connectivity port	No
Dimensions	160 x 280 x 85 mm
Power	200-240V / 50-60Hz / 20W

SciQuip GyroStir D – Digital 5” Magnetic Stirrer

- Stainless Steel 5” plate
- Speed range 100 – 1,500 rpm
- Max stirring capacity 20 litres (water)
- LCD display – both set & real time values
- Maintenance free brushless motor
- RS232 connectivity interface port
- Supplied with 1x stirring bar



Ordering Information

Part No.	Description
SP2230-D	GyroStir D - Digital Magnetic Stirrer, Stainless Steel Plate

Technical Specifications

Workplate Dimensions	Ø 135 mm
Max Stirring volume (H2O)	20L
Speed Range	100 – 1,500 rpm
Display	LCD
Connectivity port	RS232
Dimensions	160 x 280 x 85 mm
Power	200-240V / 50-60Hz / 20W

SciQuip GyroStir 10 – Analogue 10 place Magnetic Stirrer

- Up to 10 positions simultaneously
- Speed Range 0 – 1,100 rpm
- Max stirring capacity 400ml (water) per position
- Individually controlled positions for constancy across samples
- Maintenance free brushless motor
- Stainless Steel plate with silicone film – anti-slip & anti-corrosion
- Supplied with 10x stirring bar



Ordering Information

Part No.	Description
SP2240-10A	GyroStir 10 Place Magnetic Stirrer, stainless steel plate

Technical Specifications

Workplate Dimensions	180 x 450 mm
Max Stirring volume (H2O)	400 ml (per position)
Speed Range	0 – 1,100 rpm
Display	Analogue Scale
Dimensions	182 x 552 x 65 mm
Power	200-240V / 50-60Hz / 20W

SciQuip GyroStir D HOT – Digital 5” Hotplate Magnetic Stirrer

- Stainless Steel 5” hotplate
- Optional external temperature probe
- Probe auto-detection & auto display mode switch
- LCD display – both set & real time vales
- Reaches temperature rapidly
- Maintenance free brushless motor
- Max plate temperature 360°C
- Display shows ‘HOT’ until below 50°C, even when switched off.
- RS232 connectivity interface port
- Supplied with 1x stirrer bar



Ordering Information

Part No.	Description
SP2235-DH	GyroStir D - Digital Hotplate Stirrer, Stainless Steel Plate

Accessories

Part No.	Description
SP2201	Temperature sensor for digital hotplate model, 230mm, used with GyroStir D / 550H
SP2202	Temperature sensor with glass coating, 230mm, used for digital hotplate model used with GyroStir D / 550H
SP2203	Support clamp for temp sensor

Technical Specifications

Workplate Dimensions	Ø 135 mm
Max Stirring volume (H2O)	20L
Speed Range	100 – 1,500 rpm
Display	LCD
Temperature Safety Limit	360 °C
‘HOT’ Warning	> 50 °C
External Temp Sensor	SP2201 / SP2202
Connectivity port	RS232
Dimensions	160 x 280 x 85 mm
Power	200-240V / 50-60Hz / 550W

SciQuip GyroStir 550H – Digital 7” Hotplate Magnetic Stirrer

- Glass-Ceramic hotplate – highest corrosion resistance
- Optional external temperature probe
- Probe auto-detection & auto display mode
- LCD display – both set & real time values
- Reaches temperature rapidly
- Maintenance free brushless motor
- Max plate temperature 580°C
- Display shows ‘HOT’ until below 50°C, even when switched off.
- RS232 connectivity interface port
- Supplied with 1x stirrer bar



Ordering Information

Part No.	Description
SP2245-550H	GyroStir 550H - Digital 7" Square Hotplate Stirrer, Glass Ceramic Hotplate

Accessories

Part No.	Description
SP2201	Temperature sensor for digital hotplate model, 230mm, used with GyroStir D / 550H
SP2202	Temperature sensor with glass coating, 230mm, used for digital hotplate model used with GyroStir D / 550H
SP2203	Support clamp for temp sensor

Technical Specifications

Workplate Dimensions	184 x 184 mm
Max Stirring volume (H2O)	20L
Speed Range	100 – 1,500 rpm
Display	LCD
Temperature Safety Limit	580 °C
‘HOT’ Warning	> 50 °C
External Temp Sensor	SP2201 / SP2202
Connectivity port	RS232
Dimensions	215 x 360 x 112 mm
Power	200-240V / 50-60Hz / 1050W

SciQuip GyroStir 10 H – Analogue 10 place Hotplate Magnetic Stirrer

- Plate temperature up to 100°C
- 10 positions
- Uniform temperature distribution
- Individually controlled positions
- Maintenance free brushless motor
- Stainless Steel plate with silicone film – anti-slip & anti-corrosion
- Supplied with 10x stirrer bars



Ordering Information

Part No.	Description
SP2240-10AH	GyroStir 10 Place Hotplate Stirrer, stainless steel plate with silicone film

Technical Specifications

Workplate Dimensions	180 x 450 mm
Max Stirring volume (H2O)	400 ml (per position)
Speed Range	0 – 1,100 rpm
Display	Scale
Temperature Safety Limit	140 °C
'HOT' Warning	N/A
External Temp Sensor	N/A
Connectivity port	N/A
Dimensions	182 x 552 x 65 mm
Power	200-240V / 50-60Hz / 490W

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.