

myFUGETM 5



- For quick spin downs of **5ml** centrifuge tubes
- Includes adapters for 12 x 75mm tubes
- Starts & stops with closing/opening of the lid
- Conserves valuable bench space

CE



Combination rotor accepts
1.5/2.0ml & 5ml centrifuge tubes



New 5ml Centrifuge Tubes
(Item: C1005-T5)

The myFugeTM 5 microcentrifuge includes a four position rotor for compatibility with a variety of 5ml test tubes, including 5ml snap cap centrifuge tubes. In addition, this combination rotor has four additional slots for centrifuging standard microtubes (1.5/2.0ml). Simply close the lid and the centrifuge quickly accelerates to the fixed speed of 5,500rpm. Open the lid, and the rotor decelerates to a stop for removing samples.

The 5ml rotor positions accepts 5ml centrifuge tubes (snap cap) directly. While 5ml blood and culture tubes (12x75mm) and cryovials are also accepted with the included adapters. The microtube positions are designed to accept 0.5 and 0.2ml tubes with the use of optional adapters.

Technical Data:

Speed:	5,500rpm / 2,000 xg
Capacity:	4 x 5ml tubes & 4 x 1.5/2.0ml
Dimensions:	5.5 x 7.9 x 5.5 in. 14 x 20 x 14 cm
Weight:	2.2 lbs. / 1 kg
Electrical:	115V 50Hz or 230V, 60Hz
Warranty:	1 Year



Ordering Information:

C1005*	myFuge 5 Mini Centrifuge with combination rotor, 4 x 5ml & 4 x 1.5/2.0ml
C1008-A5-6	0.5ml Adapter Pack, 6/pk
C1008-A2-6	0.2ml Adapter Pack, 6/pk
C1005-AC2	Adapter Pack, for Cryovials (0.5 to 2.0ml) and 12x32mm HPLC vials, 4/pk
C1005-T5	5ml centrifuge tube, clear, 200/pk (2 bags of 100)
C1005-T5-ST	5ml <i>sterile</i> centrifuge tube, clear, 200/pk. (10 bags of 20)

*115V with US plug. To order in 230V, add (-E) to the item number.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.