



# CAPP *Controller*

## MOTORIZED PIPETTE CONTROLLER

### Features:

- Light and fatigue free: weighs only 167g including the battery
- Lithium battery with no memory effect and up to 12 hours of continuous use
  - Can be charged while in use
- Powerful: 50mL pipette is filled in less than 10 seconds
  - Express dispensing: depress the dispensing button beyond the built-in point of resistance
- Gravity delivery and blow-out dispensing mode selection
  - Fully adjustable speed-setting
- Exchangeable 0.2  $\mu$ m PTFE disc filter acts as effective shut-down system if overfill occurs
  - Exchangeable pipette holders are made of silicon
- For use with glass or plastic pipettes ranging from 0.1mL to 100mL
  - Vapour outlet for venting out corrosive liquids or fumes
- Fully autoclavable nose cone
  - UV-resistant housing
- Complies with CE and RoHS



## Ordering Information

Cat. No.	Description
PA-100	CAPPController Motorized Controller, 0.1-100mL



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.