



Display indicative only

C2006 Small Dual Centrifuge

C2006. (230V 50/60Hz). 1.C2006. (110V 60Hz)

| | |
|-------------|---------------------------------|
| Speed max | 500-6,000 Rpm (10 Rpm steps) |
| Rcf max | 3,800 G |
| Timer | 0-99 Mins & Hold (30 sec steps) |
| Dims HWD | 275 x 325 x 470mm |
| Weight | 21.8Kg (without rotor) |
| Power | 160 watts |
| Memory | 10 programs |
| Accel Rates | 10 programs |
| Decel Rates | 10 programs |

Rotors 6000 Rpm



| Rotor | BRK5412 | BRK5406 |
|--------------------------|------------|------------|
| Rotor type | 12 x 15ml | 6 x 50ml |
| Tube size max | 17 x 120mm | 30 x 120mm |
| Minimum speed Rpm | 500 | 500 |
| Maximum speed Rpm | 6,000 | 6,000 |
| Maximum rcf (G) | 3,800 | 3,800 |
| Radius max cms | 9.5 | 9.5 |
| Sample tube angle ° | 30 | 30 |
| Acceleration time (secs) | 25 | 25 |
| Deceleration time (secs) | 25 | 25 |
| Autoclavable (frequency) | 121°C (20) | 121°C (20) |

Reducers (Pack of 4)



BRK5412

BRK5406

| Rotor | BRK5412 | BRK5406 |
|---------------|-------------|-------------|
| Part number | RM05 (5ml) | RM15 (15ml) |
| Tube size max | 13 x 80mm | 17 x 120mm |
| Part number | RM10 (10ml) | RM25 (25ml) |
| Tube size max | 13 x 100mm | 25 x 100mm |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.