

## Ceti Magnum Trinocular Compound Microscope (4x, 10x, 40x, 100x Plan Achromatic)

**Product Code:** 2722.0000M

**Brand:** Ceti

Ceti's Magnum Trinocular Compound Microscope (4x, 10x, 40x, 100x Plan Achromatic) is a popular, economically priced, biological upright microscope. Providing excellent quality images via the wide-field FN20mm eyepieces, this microscope includes four Semi Plan Achromatic objectives in the reversed quintuple nosepiece.

It features a robust, T-shaped base, designed to offer maximum stability. Modern, ergonomic features such as a 30° inclined eyepiece tube, ensure a comfortable viewing position over a prolonged period.

This model can also be used in Brightfield mode so it is extremely flexible. Bright, white LED transmitted lighting, with intensity control, is also featured.

### Features

- Sturdy base for comfortable viewing over long periods
- Wide-field 10X, FN20 eyepieces with dioptic adjustment
- Quintuple, reversed nosepiece with four Plan Achromatic objectives
- Trinocular eyepiece tube as standard to accept digital microscope cameras
- Bright, white LED illumination with intensity control
- Full Köhler illumination



## Specifications

<b>Optical Head</b>	Trinocular with photo port, inclined at 30°
	Interpupillary Distance Adjustment (Siedentopf Type): 55 to 75mm
	Diotropic adjustment on one eyepiece
<b>Eyepieces</b>	10x/20mm wide-field
<b>Objectives</b>	4x, 10x, 40x, 100x Plan Achromatic
<b>Condenser</b>	Pre-condenser with field diaphragm
	Brightfield Abbe - N.A. 1.25
	Focusable and centrable, iris diaphragm
<b>Stage</b>	210 x 140mm with built-in mechanical stage and slide holder
	90 x 60mm, Vernier scale
<b>Focusing Knobs</b>	Coaxial coarse and fine with stop
	Tension adjustment ring
<b>Mains Power</b>	220 - 240V/50 - 60Hz
<b>Illumination, Transmitted</b>	3W LED with intensity control
<b>Supplied With</b>	Blue and green filter, 2 spare fuses, immersion oil (10ml), soft rubber eyecups, dust cover, and instruction manual



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.