

Ceti Versus-TC LED Trinocular Inverted Microscope with Integrated Phase Contrast

Product Code: 3660.4000M

Brand: Ceti

Ceti's Versus-TC LED Trinocular Inverted Microscope with Integrated Phase Contrast is a robust tissue culture microscope, equipped with three phase contrast objectives as well as a 10x brightfield objective that is useful for visualising a larger field of view. The phase slider is fitted with centrable phase rings plus an open position for brightfield viewing.

The Versus-TC LED also features a trinocular head with a camera port to accept a wide variety of cameras, a pair of 10x super wide-field (FN22) eyepieces, a quintuple position nosepiece with convenient click stop so you know that the selected objective is aligned correctly, infinity corrected, long working distance plan achromatic optics (10x brightfield in addition to 10x, 20x, and 40x phase contrast), a large mechanical stage equipped with a range of Petri dish and culture flask holders, ergonomic coaxial coarse and fine focus and a bright white, adjustable intensity, 3W LED light.



Features

- Compact, sturdy base for many years of use
- Super wide-field 10X, FN22 eyepieces with dioptic adjustment for a wide field of view
- Quintuple, reversed nosepiece with four LWD Plan Achromatic objectives
- Eyepiece tube with photo port to accept digital microscope cameras
- Bright, white LED illumination with intensity control
- Phase contrast condenser with centrable phase rings

Specifications

Optical Head	Trinocular with photo port, inclined at 45°
	Interpupillary distance adjustment (Siedentopf type): 48 to 75mm
	Dioptic adjustment on both eyepiece tubes
Eyepieces	10x/22mm wide-field, centring telescope
Nosepiece	Quintuple reversed, revolving with click stop
Objectives	10x Plan Achromatic Brightfield LWD 10x, 20x, 40x Phase Contrast Plan Achromatic LWD
Condenser	Phase contrast condenser with adjustable centering and brightfield position
	Working Distance: 70mm
	Focusable and centrable, iris diaphragm
Stage	300 x 170mm with additional mechanical stage
	75 x 132mm movement, Vernier scale
Focusing Knobs	Coaxial coarse and fine with stop
	Tension adjustment ring
Mains Power	220 - 240V/50 - 60Hz
Illumination	3W LED with intensity control
Supplied With	Blue and green filters, 2 spare fuses, immersion oil (10ml), Allen key, dust cover, and instruction manual



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.