

## Clifton Targa Centrifuge NE 010G

*Clifton Range®*



Targa Centrifuge is the latest model in the Clifton Range designed to suit a variety of applications in multiple markets. It is compact and robust, making it ideal for routine low speed centrifugation. Perfect for fractionation of samples and simple separation of biological cells.

Constructed with a durable steel body. The rotor is removable for easy cleaning and push button control panel is wipe clean.

### NE 010G

GENERAL	
Rotor Type	Supplied with Angle Head Rotor, 6 aluminium buckets and 6 x 15ml clear tubes.
Capacity	Rotor suitable for tube capacity: 6 x 15ml maximum inclined 45°. Rotor tube dimensions: 18mm diameter bucket, maximum tube length 125mm or 120mm conical capped. Accessory tube adaptors for 3ml, 4ml, 5ml tubes
Maximum rotor speed	4000rpm
Maximum RCF	2075g
Maximum radius	130mm
Rotor acceleration	<20 secs max speed
Rotor imbalance cut off	Yes
Lid lock safety	Yes
Lid drop protection	Yes
Lid release	Auto lid release with audible sounder when rotor stops
Noise level	<52dB
Timer	Yes
Program memory	10
Operation	Digital LED display with rotor speed in rpm or rcf and timer in minutes or seconds. Start stop, program and lid release buttons.
Power	440W
Voltage	230V 50/60Hz
DIMENSIONS	
External (mm)	296w x 412d x 206h mm
Packed Weight	19 kg

# Clifton Targa Centrifuge

## Accessories



### Tube Adaptor 3ml (Rotor accepts 6 adaptors)

Accessory Code	Description
T010-3	With adaptor added to rotor buckets 3ml tubes can be centrifuged



### Tube Adaptor 4ml (Rotor accepts 6 adaptors)

Accessory Code	Description
T010-4	With adaptor added to rotor buckets 4ml tubes can be centrifuged



### Tube Adaptor 5ml (Rotor accepts 6 adaptors)

Accessory Code	Description
T010-5	With adaptor added to rotor buckets 5ml tubes can be centrifuged




### Tube 15ml Conical Tube

Accessory Code	Description
429910	15ml Conical tube, unskirted and capped (Pack of 20)



### Tube 15ml Conical Tube

Accessory Code	Description
429946	15ml Conical tube, unskirted, capped, sterile and bagged (Pack of 20)



### Cap - Push on, for 16mm Tubes

Accessory Code	Description
307202	Internal ribs to clamp firmly on tube, leakproof (Pack of 20)



### Round Bottom Tube (Packs of 20)

Accessory Code	Description
	400900 : Clear polypropylene 16mm x 95mm, graduated 2.5, 5 and 10 ml
	400705 : Clear polypropylene 16mm x 100mm
	300705 : Clear polystyrene 16mm x 100mm
	300900 : Clear polystyrene 16mm x 95mm, graduated 2.5, 5 and 10ml



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.