



Model: CMCF650

650L Laboratory Chest Freezer

Overview

- -10°C to -25°C Temperature Range
- External Digital Control and Temperature Display
- Fitted Handle with Key Lock and 2 Keys
- Twin Baskets to Separate Frozen Items
- 4 Star Freezer Rating

5 YEARS PARTS & LABOUR WARRANTY

ANNUAL SERVICE & CALIBRATION AVAILABLE

Key Features



Controller

Controller system allows reading, recording, and resetting with the touch of one button.



Flexible Storage

Multiple baskets to provide a flexible storage solution.



Spark Free

No internal light or fan, suitable for a spark free laboratory.

Details

Model: CMCF650

Gross Capacity (litres)	567	Temperature Display on Panel, Temp Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms	Yes
Net Capacity (litres)	560	Number of Shelves/Baskets	2
Unit External Width (mm)	1805	Defrost Drain	With
Unit External Depth (mm)	750	Shelf Material	White Wire
Unit External Height (mm)	850	Interior Light	No
Product Weight (kg)	79	Light Type	N/A
Packed Width (mm)	1850	Fan	No
Packed Depth (mm)	780	Spark Free	Yes
Packed Height (mm)	880	Refrigerant Type	R600a
Packed Weight (kg)	86	CFC Free, HCFC Free, Ammonia Free	Yes
Lid Thickness (mm)	80	Voltage Supply	220-240V/ 50Hz
Wall Thickness (mm)	70	Certification	CE/UK
Temperature Range	-10°C to -25°C	Input Power (W)	125
External Finish, Internal Finish	Aluminium, PCM	Energy Consumption (KWh/annum)	314
Battery Backup	No	Energy Rating	E
Door Type, Door Style	Flat, PCM, Hinged		
Open Door Sensor	No		
Lock	Yes, Integrated with handle		
Adjustable Feet (number)	4		
Cooling System	Static		

Related Products



CMF86V208
-86°C Ultra Low
Temperature Freezer



CMF86V80
-86°C Ultra Low
Temperature Freezer



CMCF550
550L Laboratory
Chest Freezer



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.