



# Compressor Station

## 091-12 (GBWA038ABK)

- Oil free compressor technology
- Maintenance free except for the intake filter
- Ready for connection
- Adjustable regulator on main control panel
- Can be used as a wall mounted or base mounted unit
- Low noise level 55 – 56 dB (A)
- Suitable for use in an S1 condition up to 30°C ambient
- Extended warranty available (on request)
- Auxiliary connection to enable external 10ltr pressure vessel to be added if required



### Technical data

Delivery Rate @ 0 bar	35 l/min
Delivery Rate @ 5 bar	20 l/min
Receiver volume	10 litres
Voltage	230 V
Frequency	50/60 Hz
Protection class	IP54
Dimensions in cm (H x W x D)	41x60x25
Weight	19.5 kg (Unpacked)



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.