

Silenced Compressor Station

280-18 (GBJA062ABKM)



- Oil free compressor technology
- Maintenance free except for the intake filter(s) and sinter filter
- Ready for connection
- Adjustable regulator on main control panel
- Internal membrane drier suitable for air purity up to (1.3.1)
- Contained within a temperature controlled acoustic cabinet
- Low noise level 59-60 dB (A)
- Suitable for use in an S1 condition up to 30°C ambient
- Extended warranty available (on request)
- Auxiliary connection to enable external pressure vessel to be added if required
- Trolley mounted for ease of transport. The ultimate portable compressor.



Technical data

Delivery Rate @ 0 bar	40 l/min
Delivery Rate @ 5 bar	28 l/min
Receiver volume	0.25 litre
Voltage	230 V
Frequency	50/60 Hz
Protection class	IP20
Dimensions in cm (H x W x D)	63 x 42 x 51
Weight	25 kg (Unpacked)



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.