



Logiclave LABC Range of Autoclaves

A range of freestanding autoclaves available in standard cylindrical chamber sizes from 100 to 200 litres capacity and intended to suit a wide variety of modern laboratory applications. This range compliments our range of rectangular section autoclaves and utilises the same proven and easy to use Logiclave type control system.

- Mounted on swivel castors for easy transportation and positioning
- Designed to fit through standard width doorways
- Hygienic stainless steel fascia panels throughout entire range

Typical Applications:

- Porous / Fabrics
- Equipment & Utensils
- Empty Glassware
- Culture Media
- Plastic Discard
- Fluid Discard

Standard Features:

- Logiclave Controller with Touchscreen Interface
- Front Loading (Waist height)
- Power Operated Vertical Sliding Door
- Floating Door Seal
- 316L Stainless Steel Chambers
- Chamber Pressure Gauge
- Hygienic Stainless Steel Enclosure
- Emergency Stop / Door Safety Interlocks

Nominal Chamber Sizes(Cylindrical Vessels):

| Litres | Diameter | Depth |
|--------|----------|-------|
| 100 | 500 | 500 |
| 145 | 500 | 750 |
| 200 | 500 | 1000 |

Overall Nominal Dimensions:

| Litres | Height | Width | Depth |
|--------|--------|-------|-------|
| 100 | 1735 | 730 | 1000 |
| 145 | 1735 | 730 | 1100 |
| 200 | 1735 | 730 | 1350 |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.