

# Infrared Heater IRB 1



## Description

The compact laboratory heater IRB 1 with a heating area of 60 x 60 mm is optimally suited for fast heating tasks without open flame up to max. 700°C.

### Basic equipment

Basic device with support cross

## Details

<b>Max. temperature:</b>	700°C
<b>Heating power:</b>	250 W
<b>Heating area:</b>	60 x 60 mm
<b>Max. ambient temperature:</b>	30°C

**Order no.:** 6060 000  
**Accessory:** Power control LR 1  
Support rod for stand connection

## **Technical specifications**

**Electrical supply:** 230 V, 50 / 60 Hz (115 V version on request)

**Dimensions (w x d x h):** 100 x 105 x 105 mm

**Weight:** 0,7 kg



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.