

# Infrared Heater IRB 6



## Description

With the series heater IRB 6 the samples can be heated on 6 separate heating areas of 60 x 60 mm each up to 700°C without open flame. Each heater can be individually used and controlled.

### Basic equipment

6 x IRB 1 with support crosses, 1 x power control unit LR 6

## Details

<b>Max. temperature:</b>	each 700°C
<b>Heating power:</b>	6 x 200 W
<b>Heating area:</b>	6 x 60 x 60 mm
<b>Power control:</b>	6 stepless controls between 0 - 100 %

**Max. ambient temperature:** 30°C

**Order no.:** 9066 000

**Accessory:** Support base  
Horizontal PVC rods for support base  
Clamps D50 (NS 45)  
Clamps D32 (NS 29)  
Clamps D20 (NS 19)  
Clamps D25 (NS 15)

## Technical specifications

**Electrical supply:** 230 V, 50 / 60 Hz (115 V version on request)

**Heating power:** max. 1200 W

**Dimensions (w x d x h):** 697 x 305 x 122 mm

**Weight:** 9 kg



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.