


PDS No.7101xx / 7101xx-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 11	Pipettes, CELLSTAR®, 2 ml				 greiner bio-one
	Greiner Item-No. 7101xx / 7101xx-TRI				
Valid for Item-No.:	710107 (sterile)	710160 (sterile)	710180 (sterile)	710183 (sterile)	
		710160-TRI (sterile)		710183-TRI (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, 2 ml, sterile 710107, -160, -180: <u>Serological Pipettes</u> : graduation (scale 1/100 ml), lateral green stripe, synthetic plug 710160-TRI: <u>Serological Pipettes</u> : graduation (scale 1/100 ml), lateral green stripe, synthetic plug, triple packed 710183: <u>Aspiration Pipette</u> : standard (without plug, graduation, lateral stripe) 710183-TRI: <u>Aspiration Pipette</u> : standard (without plug, graduation, lateral stripe), triple packed
1.2	Dimensions	See Customer Drawing
1.3	Volume	2 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	Pipette: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Plug: white
1.6	Sterilization	710107, -160, -180, -183: SAL 10 <sup>-3</sup> 710160-TRI, 710183-TRI: SAL 10 <sup>-6</sup>
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	710107	710160	710180, -183	710160-TRI,	710183-TRI
3.1	Pieces / Bag	25 Bulk packs (plastic)	1 / 100 "peel off"-packaging (plastic/plastic)	1 / 100 "peel off"-packaging (paper/plastic)	1 "peel off"-packaging (plastic/plastic)  10 "Triple packed"	1 "peel off"-packaging (paper/plastic)  10 "Triple packed"
3.2	Pieces / Box	1000	1000	1000	100	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)				
3.4	Other Information	Certificate of Quality				

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 10	Date 1 February 2018	Date 2 February 2018	Date 2 February 2018	
Date 14.08.2017	Name S. Kaelberer	Name Dr. J. Bischoff	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.