

## Pharmacy Refrigerators HYC-360



### Five Alarm Functions

High/low temperature, power failure, door ajar, sensor error, low battery



### Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



### Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements

### Ergonomic Design



High Accuracy Temperature Control



Interior Fluorescent Light



Lockable Design

# Specifications

## HYC-360

### Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	55

### Performance

Temp Range(℃)	2~8
---------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220~240/50	208~230/60	115/60
Power(W)	315	315	340
Electrical Current(A)	2.2	2.2	4.5

### Dimensions

Capacity(L/Cu.Ft)	360/12.7
Net/Gross Weight(approx)	99/114(kg) 218.3/251.3(lbs)
Interior Dimensions(W*D*H)	550*460*1340(mm) 21.7*18.1*52.8(in)
Exterior Dimensions(W*D*H)	620*655*1995(mm) 24.4*25.8*78.5(in)
Packing Dimensions(W*D*H)	695*720*2080(mm) 27.4*28.3*81.9(in)
Container Load(20'/40'/40'H)	24/48/48



### Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### Accessories

Caster	Y
Foot	Y
Shelves /Basket	5/-

### Others

Certificate	CE	/
-------------	----	---



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.