



## Rapid, accurate colony counter

The ecount™ is a simple to use colony counter for counting cells in microbiology applications or tissue culture. This unit combines a multi-function electronic counter with a Sharpie® pen for counting and marking to prevent missing or double-counting colonies in a Petri dish.

- US Patent No. D564385 and EP Registered Design No. 705751-0001 featuring a light, ergonomic body that is balanced for easy manipulation
- The unit has convenient easy to use touch controls for multiple counts, sums, counting up or down and turning the unit on/off
- Count up/down verification, enables you to confirm you counts
- Easily verify each count with either an audible beep or LED flash
- Displays totals for up to 32 individual counts or zones
- Quick and easy to change between different colored markers to check different colony counts in the same view
- Multi-functional counter can be turned off allowing the unit be used for writing
- Up to 999 colonies per zone helps simplify the counting process
- Long Life – CR1225 3.0 V Lithium battery

Part No. **HS120000**

### Specifications

**Dimensions:** Unit is 3 cm (1.2 in) wide and 19 cm (7.4 in) long (with pen)

**One-year warranty**

Sharpie® is a registered trademark of Sanford, A Newell Rubbermaid Company



Weighted stand for bench top use



Thumb button allows you to use the ecount™ as a digital tally counter



Heathrow Scientific®



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.