



5
Year
Warranty

Back-lit Display
showing time, speed and
operational information on
the same screen

Adjustable
time and speed

Quick-spin Button
for momentary spin
requirements, press and hold to
operate, release to end spin cycle

EASY SWITCH OVER
NO
TOOLS REQUIRED

Gusto® High-Speed Mini Centrifuge

Maximum speed with minimal space

This compact centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement

Made of durable ABS with UV-Resistant coating

User Friendly

- Intuitive Push button programming and bright LCD screen make it easy to use

Quiet Operation

- Brush-less motor enables quiet operation and low vibration

Safety Features

- Lid safety feature prevents opening while rotor is turning and clear polycarbonate allows user to see rotor movement for added visual safety
- Includes automatic shut-off feature for sample imbalance and tilt

Specifications

Speed Range	500–12,500 rpm
Run Time Ranges	1 to 30 minutes
RCF	15 - 9,783 x g (max RCF calculated)



Package Includes

- Gusto® High-Speed Mini Centrifuge
- 12-Place 1.5/2.0 mL Tube Rotor
- 32-Place 0.2 mL Tube Rotor (singles, or 4 x 8 Strip)
- 12 Tube Adapters 0.2 mL
- 12 Tube Adapters 0.5 mL
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Item No.	Color	L x W x H in	L x W x H cm
HS10050	Blue ■	8 x 6.7 x 4.5	20.3 x 17.1 x 11.4

Replacement parts

HS100502	Universal Power Cord Set
HS100503	Gusto® Rotor 12 x 1.5/2.0 mL Tube (with Cover and Knob)
HS100504	Gusto® Rotor 4 x 8 Strip 0.2 mL Tube (with Knob)
120015	Tube Adapters, 0.2, 0.5 mL, Pack of 12 (6 ea.)
HS100506	Gusto® Rotor Cover (for 12 x 1.5/2.0 mL Tube Rotor Only)
120527	Power Adapter



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.