

SE Range

Low Temperature Chest Freezer

TEFCOLD 

New product

Features and Benefits

- * Electronically Controlled
- * Lock Fitted
- * Low Temperature -45°C
- * No Interior Light
- * 80mm Insulation
- * Defrost Drain

Technical Specifications

Temperature Range . . . -40/-45°C
Refrigerant R404A
Type of Defrost Manual
No. of Baskets 1
Exterior Finish White
Interior Finish Aluminium
Power Required 13 Amp Supply

SE40



SE30

Dimensions	Height:	Width:	Depth:	Weight:	Capacity:
External					
SE10	895mm	734mm	705mm	46kg	152ltr. (5.37cu.ft.)
SE20	895mm	984mm	705mm	51kg	230ltr. (8.12cu.ft.)
SE30	895mm	1284mm	705mm	61kg	323ltr. (11.4cu.ft.)
SE40	895mm	1504mm	705mm	71kg	392ltr. (13.84cu.ft.)
Internal					
SE10	685mm	570mm	440mm	2	Max No of Baskets
SE20	685mm	820mm	440mm	3	
SE30	685mm	1120mm	440mm	5	
SE40	685mm	1340mm	440mm	6	

MODEL	Description	MODEL	Description
SE10	White	SE30	White
SE20	White	SE40	White

LT Range

Low Temperature Chest Freezer


Elcold 

Features and Benefits

- * Electronically Controlled
- * Lock Fitted
- * Super Low Temperature to -45°C
- * Inner Liner of Hammered Aluminium
- * ISO9001 Quality Controlled Manufacture
- * 100mm Insulation
- * Digital Temperature Display
- * Temperature Alarm (Visual Only)
- * No Interior Light

Technical Specifications

Temperature Range . . . -25°C/-45°C
Refrigerant R404A
Type of Defrost Manual
No. of Baskets supplied
11LT 1
21LT 1
31LT 1
41LT 3
51LT 3
Exterior Finish White
Interior Finish Aluminium
Power Required 13 Amp Supply



Dimensions	Height:	Width:	Depth:	Weight:	Capacity:
External					
EL11LT	860mm	720mm	730mm	46kg	136ltr. (4.8cu.ft.)
EL21LT	860mm	1050mm	730mm	58kg	237ltr. (8.4cu.ft.)
EL31LT	860mm	1300mm	730mm	68kg	314ltr. (11.1cu.ft.)
EL41LT	860mm	1500mm	730mm	82kg	375ltr. (13.3cu.ft.)
EL51LT	860mm	1700mm	730mm	90kg	438ltr. (15.47cu.ft.)
Internal					
EL11LT	660mm	520mm	450mm	1	Max No of Baskets
EL21LT	660mm	850mm	450mm	2	
EL31LT	660mm	1100mm	450mm	3	
EL41LT	660mm	1300mm	450mm	4	
EL51LT	660mm	1504mm	450mm	4	

MODEL	Description	MODEL	Description
EL11LT	White	EL41LT	White
EL21LT	White	EL51LT	White
EL31LT	White		

FROZEN STORAGE



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.