



## Safety



Fire simulation test (90mins)



Fire performance test (NFPA 30 Vs. EN 14470-1)

### Fire Safety Storage Cabinet

- Safe storage of flammable / combustible materials in the event of a surrounding fire.
- Designed according to the strictest European standards EN 14470-1
- Fire proof performance : 30 min / 90 min

Internal Volume	Tall type (90min type) : 763L / 357L Tall type (30min type) : 802L / 375L Underbench type (90min type) : 244L / 110L Underbench type (30min type) : 252L / 116L
Model	SC3-F + EM-33

※ Available to install Safety Waste Disposal System. (Option, p43)

## Safety



### Acid Corrosive Cabinet

- Acid/Base reagents can be safely stored for a long period of time.
- Solid steel structure with corrosion-resistant interior finish.

Internal Volume	Stand Type : 894L / 852L / 445L / 426L Under Type : 314L / 227L / 146L / 115L
Model	SC-C + EM-33

※ Available to install Safety Waste Disposal System. (Option, p43)

### Molded Acid Corrosive Cabinet

- Durable chemical resistant plastic storage cabinets for Acid/Base materials reagents.
- Easy use with connecting filter or fan modules.

122L  
Exterior Dimension : 720 x 472 x 535 mm

PCC-720

## General



### Multipurpose Cabinet Window

- Window door type for easy figures out stored samples.
- Drawer shelves with PP having good chemical -resistance, makes easy store.
- Fume extraction module & filters. (option)

Model DS-OW

## Functional



### Filtering Storage Cabinet

- Filter installed inside, internal air circulation.
- Easy moving and installation, safety warning alarm.

• Internal Volume: 604 L / 1317 L  
• Capacity (L or 1L Bottle): 140 / 280  
• Optional Filter : HEPA Filter, Complex Carbon Filter (Aldehyde, acid and amines are available.)

FSC



### Gas Cylinder Cabinet

- Safe protection from external shock.
- Convenient folding ramp structure for replacing heavy cylinders.

Gas Cylinder Max. storage capacity (Ø230 mm) : 2 ea / 3 ea

DS-GC



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.