


OFI 202-40M (oil-less)

Voltage	V	230	230	3x400	3x400								
Frequency	Hz	50	60	50	60								
Motor	HP	2.0	2.0	2.0	2.0								
	kW	1.47	1.47	1.47	1.47								
Displacement	l/min	290	328	290	328								
	CFM	10.24	11.58	10.24	11.58								
FAD @ 8 bar/120 psi	l/min	130	146	130	146								
	CFM	4.59	5.16	4.59	5.16								
Max pressure ¹⁾	bar	8	8	8	8								
	psi	120	120	120	120								
Max. current	A	8	9	5	5								
Tank size	litres	40	40	40	40								
	gallon	10.6	10.6	10.6	10.6								
Weight	kg	116	116	116	116								
	lbs	256	256	256	256								
Noise level	db(A)	60	63	60	63								
Thermal protection		Yes	Yes	Yes	Yes								
Duty cycle		100%	100%	100%	100%								



Dimensions (l x w x h)
661 x 643 x 856 mm / 26 x 25 3/8 x 33 6/8 inch

¹⁾ Higher pressure available upon request

Technical modifications reserved

JUN-AIR 

IDEX
IDEX CORPORATION





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.