


2xOFI 202-90B (oil-less)

Voltage	V	200	200	230	240	3x400	3x440					
Frequency	Hz	50	60	50	60	50	60					
Motor	HP	4.00	4.00	4.00	4.00	4.00	4.00					
	kW	2.94	2.94	2.94	2.94	2.94	2.94					
Displacement	l/min	580	656	580	656	580	656					
	CFM	20.48	23.17	20.48	23.17	20.48	23.17					
FAD @ 8 bar ³⁾	l/min	260	292	260	292	260	292					
	CFM	9.18	10.31	9.18	10.31	9.18	10.31					
Max. pressure ¹⁾	bar	8	8	8	8	8	8					
	psi	120	120	120	120	120	120					
Max. current	A	-	-	16.0	18.0	10.0	10.0					
Tank size	litres	90	90	90	90	90	90					
	gallon	23,7	23,7	23,7	23,7	23,7	23,7					
Weight	kg	115	115	115	115	115	115					
	lbs	253	253	253	253	253	253					
Noise level	dB(A)/1m	79	78	79	78	79	78					
Thermal protection		Yes	Yes	Yes	Yes	Yes	Yes					
Duty cycle		100%	100%	100%	100%	100%	100%					



Dimensions (l x w x h)
1000 x 530 x 850 mm / 39.4 x 19.7 x 33.5 inch

¹⁾ Higher pressure available upon request

³⁾ Displacement is reduced by approx. 18-20% on units with dryer (D). Min. pressure required to operate dryer: 6 bar.

Technical modifications reserved

JUN-AIR 





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.