



PG102C Countertop Pharmacy Refrigerator Glass Door 41L

PG102C

- ✓ 2°C to 8°C set temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 2 year warranty as standard

SPECIFICATION

| | |
|------------------------------|--|
| Dimensions | 379mm W x 702mm H x 500mm D |
| Weight | 25.5kg |
| Product name | Pharmacy Fridge Glass Door 41L PE102C |
| GDPA code | PE102C |
| SKU code | UK 444441791, EU 444442331, AU 444442332 |
| External dimensions | 379mm W x 702mm H x 500mm D |
| Internal dimensions | 245mm W x 400mm H x 310mm D |
| Packed dimensions | 450mm W x 740mm H x 560mm D |
| Product weight | 25.5kg |
| Packed weight | 27.5kg |
| Gross capacity | 41L |
| Net capacity | 35L |
| External finish | Painted steel |
| Colour | White/green/glass |
| Internal finish | HIPS* |
| Type of plastic | SRL800 (HIPS/ABS mix) |
| Door type | Glass |
| Door finish | Painted steel |
| Door style | Hinged right hand only, non reversible |
| Door open depth | 810mm |
| Replacement door seal | N |
| Open door sensor | Y |
| Lock | Y |
| Adjustable feet | N |
| Temperature range | 2°C to 8°C |
| Temperature display on panel | Y |
| Temperature alarm | Y |
| Number of shelves | 2 |
| Shelf material | Wire |
| Number of shelf positions | 5 |

| | |
|-------------------------|---------------------------|
| Shelf size | 270mm W x 6mm H x 308mm D |
| Interior light | N |
| Auto or manual defrost | Auto |
| Fan | Y |
| Controller type | Digital, external |
| Climate class | SN |
| Refrigerant type | R600a |
| CFC free | Y |
| HFC free | N |
| HCFC free | Y |
| Ammonia free | Y |
| Optional plug types | Y |
| Voltage supply | 230 |
| Spark free | N |
| Auxiliary output | Y |
| Insulation thickness | 40mm |
| Insulation type | Cyclopentane |
| Noise level | <45dB |
| Energy consumption | 1.68 kw/24h |
| Door open alarm | Y |
| Air temperature alarm | Y |
| Antimicrobial handle | Y |
| Optional drawer system | N |
| Stacking kit compatible | N |





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.