



FREESTANDING LABORATORY FRIDGE SOLID DOOR 475L WHITE



- ✓ **+2°C to +10°C Operating Range**
Conforms to MHRA guidelines
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Castors**
For ease of movement
- ✓ **Wire Shelves**
Enables optimum air flow inside unit



SKU CODE: PRODUCT CODE:

444441801 LR1607C

DESCRIPTION:

Freestanding Laboratory Fridge Solid Door 475L White



FULL SPECIFICATION – LR1607C FREESTANDING LABORATORY FRIDGE SOLID DOOR 475L WHITE



| | |
|----------------------------------|-----------------------|
| Description | Freestanding |
| SKU Code | 444441801 |
| EAN Code | 5034648485696 |
| GDPA Code | LR1607C |
| Gross Capacity (L) | 475 |
| Net Capacity (L) | 441 |
| Unit External Width (mm) | 695 |
| Unit External Depth (mm) | 910 |
| Unit External Height (mm) | 1688 |
| Product Weight (kg) | 95 |
| Packed Width (mm) | 785 |
| Packed Depth (mm) | 995 |
| Packed Height (mm) | 1865 |
| Packed Weight (kg) | 116 |
| Country of Origin | UK |
| Internal Dimensions (mm) (WxDxH) | 550 x 670 x 1170 |
| External Finish | White |
| Internal Finish | White |
| Type of Plastic | SRL800 (HIPS/ABS Mix) |
| Temperature Range | +2°C to +10°C |
| Colour | White/Green |
| Door Type | Solid |
| Door Style | Hinged, R/H only |
| Door Open Depth (mm) | 1596 |
| Replacement Door Seal | No |
| Open Door Sensor | No |
| Lock | Yes |
| Adjustable Feet | Castors |
| Temperature Display on Panel | Yes |
| Temperature Alarm | Yes |
| Number of Shelves | 5 |
| Number of Shelf Positions | 30 |
| Shelf Size | 5: 536 x 609 x 6 |
| Shelf Material | Wire |
| Moulded Drawers | 0 |
| Adjustable Door Bins | 0 |
| Interior Light | No |
| Fan | No |
| Spark Free | Yes |
| Auto or Manual defrost | Automatic |
| Controller type | Digital, External |
| Climate Class | N |
| Refrigerant Type | R134a |
| CFC Free | Yes |
| HFC Free | No |
| HCFC Free | Yes |
| UK 3 Pin Plug | Yes |
| Other Plug Types Available | Yes |
| Voltage Supply | 230 |
| Aux Output | No |
| Insulation Thickness (mm) | 65 |
| Insulation Type | Cyclopentane |
| Anti Microbial Handle | No |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.