



RLPR0217- ATEX Bench Top Sparkfree ATEX Refrigerator

| Part Number | RLPR0217-ATEX |
|-----------------------------|---|
| Capacity (Litres) | 66 |
| Style | Bench Top |
| Exterior Dimensions (HxWxD) | 740 x 450 x 510mm |
| Interior Dimensions (HxWxD) | 590 x 333 x 300mm |
| Weight | 34kg |
| Energy Consumption | 0.6kWh/24hr |
| Operating Temperature | +0°C to +10°C |
| Factory Set Point | +5°C |
| Interior Construction | Painted Aluminium |
| Exterior Construction | Powder Coated Steel |
| Door Style | Self Closing |
| Door Hinging | Right (can be reversed please ask at time of order) |
| Door Lock | ✓ |
| Door Earth Connection | ✓ |
| Condensate Collection | ✓ |
| Shelves | 3 |
| Shelf dimension (WxD) | 2 x 346 x 203mm 1 x 346 x 150mm |
| Spare shelves available | ✓ |
| Feet | 4 |
| Temperature Display | ✓ |
| High Temperature Alarm | ✓ |
| Low Temperature Alarm | ✓ |
| Door Open Alarm | ✓ |
| Off cycle auto Defrost | ✓ |
| Low GWP Refrigerant | ✓ |
| 2 Year Warranty (UK only) | ✓ |

The Labcold ATEX refrigerator range has been designed with an interior which is licensed as a Zone 2 hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

Warning: Do not store open containers of volatile substances in this fridge.

CLAIM
YOUR
FREE
2
YEAR
WARRANTY



All Labcold sparkfree ATEX refrigerators are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.