

[Click Below](#)

## Rotators

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E

**DUAB**

Your Accessible Solution

### For General Use

Manual pipettes  
Pip. Filler/ Stepper  
Electronic Dispenser  
Magnetic Stirrers  
Shakers

### For Analytical Labs

Electronic Titrator  
Manual Dispenser  
Rotary Evaporator  
Hotplate  
Overhead Stirrers

### For Life Science Labs

Electronic Pipettes  
SafeVac  
Centrifuges-LS  
Mixers / Vortex  
Thermo Mixers

### For Clinical Labs

Centrifuges-Clinical  
Rollers and Rotators  
Thermo Cyclers  
Spectrophotometers



MX-RD-Pro

MX-RD-E

MX-RL-Pro

MX-RL-E

Advantages for the user	Features
<b>Different models</b>	The available 4 types of rotators are round and linear, digital and economical
<b>One machine double uses</b>	The main body of MX-RD-Pro and MX-RL-Pro can be exchanged each other
<b>Wide application</b>	A wide variety of accessories are available

	MX-RD-E	MX-RD-Pro	MX-RL-E	MX-RL-Pro
<b>Arguments</b>				
<b>Performance W/rpm</b>	Two versions: for 10-70 (or 0-80 rpm s. below)			
<b>Tilt Angle</b>	Adjustable mixing angle between 0—90°		Mixing angle adjustment is not possible	
<b>Comfort</b>	Many accessories, up to 80 rpm	Speed control and optional timer, many accessories, LCD	Many accessories, up to 80 rpm	Speed control and optional timer, many accessories, LCD

[Click Below](#)

**For General Use**

Manual pipettes

Pip. Filler/ Stepper

Electronic Dispenser

Magnetic Stirrers

Shakers

**For Analytical Labs**

Electronic Titrator

Manual Dispenser

Rotary Evaporator

Hotplate

Overhead Stirrers

**For Life Science Labs**

Electronic Pipettes

SafeVac

Centrifuges-LS

Mixers / Vortex

Thermo Mixers

**For Clinical Labs**

Centrifuges-Clinical

Rollers and Rotators

Thermo Cyclers

Spectrophotometers

**Rotators**

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E



Your Accessible Solution

**MX-RD-Pro**

**Technical data**

Motor type	DC motor
Tilt angle	0-90°
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage	100-240V, 50/60Hz
Power	40W
Dimensions [D × W × H ]	280 × 210 × 300mm
Permissible ambient temperature and humidity	5~40°C, 80%
Protection class	IP21



**MX-RD-E**

**Technical data**

Motor type	DC motor
Tilt angle	0-90°
Speed range	0-80rpm
Speed display	Scale
Timer display	Scale
Operation mode	Continuous
Voltage	100-240V, 50/60Hz
Power	20W
Dimensions [D×W×H ]	300 × 220 × 310mm
Permissible ambient temperature and humidity	5~40°C, 80%
Protection class	IP21



[Click Below](#)

### For General Use

- Manual pipettes
- Pip. Filler/ Stepper
- Electronic Dispenser
- Magnetic Stirrers
- Shakers

### For Analytical Labs

- Electronic Titrator
- Manual Dispenser
- Rotary Evaporator
- Hotplate
- Overhead Stirrers

### For Life Science Labs

- Electronic Pipettes
- SafeVac
- Centrifuges-LS
- Mixers / Vortex
- Thermo Mixers

### For Clinical Labs

- Centrifuges-Clinical
- Rollers and Rotators
- Thermo Cyclers
- Spectrophotometers

## Rotators

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E



Your Accessible Solution

### Varieties of accessories



#### Disk Accessories Type-1

Used for MX-RD-Pro for 1.5mL x 60 centrifuge tubes holder

Cat. No.  
18900160



#### Disk Accessories Type-2

Used for MX-RD-Pro for 15mL x 16 centrifuge tubes holder

Cat. No.  
18900161



#### Disk Accessories Type-3

Used for MX-RD-Pro for 50mL x 8 centrifuge tubes holder

Cat. No.  
18900162



#### Disk support rods

4 pcs in a group use with disk accessories

Cat. No.  
18900140



#### Rotisserie Accessory Type-1

Used for MX-RL-Pro, for 1.5mL x 48 centrifuge tubes held horizontally

Cat. No.  
18900142



#### Rotisserie Accessory Type-2

Used for MX-RL-Pro, for 15mL x 24 centrifuge tubes held horizontally

Cat. No.  
18900143



#### Rotisserie Accessory Type-3

Used for MX-RL-Pro, for 50mL x 24 centrifuge tubes held horizontally

Cat. No.  
18900144



#### Rotisserie Accessory Type-4

Used for MX-RL-Pro, for 1.5mL x 32 centrifuge tubes held vertically

Cat. No.  
18900145



#### Rotisserie Accessory Type-5

Used for MX-RL-Pro, for 15mL x 16 centrifuge tubes held vertically

Cat. No.  
18900146



#### Rotisserie Accessory Type-6

Used for MX-RL-Pro, for 50mL x 16 centrifuge tubes held vertically

Cat. No.  
18900147



#### Plastic Clips

Used with MX-RL-E/MX-RD-E Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL centrifuge tubes



[Click Below](#)

**For General Use**

Manual pipettes  
 Pip. Filler/ Stepper  
 Electronic Dispenser  
 Magnetic Stirrers  
 Shakers

**For Analytical Labs**

Electronic Titrator  
 Manual Dispenser  
 Rotary Evaporator  
 Hotplate  
 Overhead Stirrers

**For Life Science Labs**

Electronic Pipettes  
 SafeVac  
 Centrifuges-LS  
 Mixers / Vortex  
 Thermo Mixers

**For Clinical Labs**

Centrifuges-Clinical  
 Rollers and Rotators  
 Thermo Cyclers  
 Spectrophotometers

# Rotators and Rollers

## Ordering Information



Your Accessible Solution

Cat. No.	P/N	Descriptions
8031312000	MX-T6-S	Classic Tube Roller, variable speed, 6 rollers
8031312100	MX-T6-Pro	LCD Digital Tube Roller, variable speed, 6 rollers
8031412100	MX-RD-Pro	LCD Digital Rotator(including a 50ml disk accessory - 18900162), used with disk accessories,variable speed, adjustable mixing angle
8031402100	MX-RL-Pro	LCD Digital Rotator (including a 50ml rotisserie accessory - 18900147),used with rotisserie accessories,variable speed, adjustable mixing angle
8031412000	MX-RD-E	Rotator(including accessory-18900327/8/9), used with disk accessories,variable speed,
8031402000	MX-RL-E	Rotator(including accessory-18900327/8/9), used with rotisserie accessories,variable speed,
<b>Accessories of Rotator</b>		
Cat. No.	P/N	Descriptions
18900160		Disk accessory, for 1.5ml x 60 centrifuge tubes holder, use with MX-RD-Pro
18900161		Disk accessory, 15ml x 16 centrifuge tubes holder, use with MX-RD-Pro
18900162		Disk accessory, 50ml x 8 centrifuge tubes holder, use with MX-RD-Pro
18900140		Disk support rods, 4 pcs, use with disk accessories
18900142		Rotisserie accessory, for 1.5ml x 48 centrifuge tubes held horizontally, use with MX-RL-Pro
18900143		Rotisserie accessory, for 15ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro
18900144		Rotisserie accessory, for 50ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro
18900145		Rotisserie accessory, for 1.5ml x 32 centrifuge tubes held vertically, use with MX-RL-Pro
18900146		Rotisserie accessory, for 15ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro
18900147		Rotisserie accessory, for 50ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro
18200995		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18200996		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18200997		Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18900327		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC
18900328		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC
18900329		Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.