



Advanced Humidity and Thermistor Data Logger

ATMDL-HUSB

- ✓ -35 to +80°C (-31 to +176°F) and 0 to 100%RH measurement range
- ✓ Stores over 16,000 readings for both temperature and humidity
- ✓ Logging rates between 10 seconds and 12 hours
- ✓ Immediate and delayed logging start
- ✓ User-programmable alarm thresholds for both temperature and humidity
- ✓ Status indication via red/green LEDs
- ✓ Environmental protection to IP67

SPECIFICATION

Dimensions	108mm W x 25mm H x 22mm D
Weight	N/A
Measurement range	-35°C to 80°C Relative Humidity = 0 to 100%RH
Temperature display on panel	Y
Accuracy	+/-0.5°C



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.