



## Upright Laboratory Control Plus Freezer 232L ATEX

GDPA code: **LSFSF232UK-ATEX ATEX UP Frz**

SKU code: **444410673**

- ✓ -18°C to -25°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in  
Medical Refrigeration

## Specification: Upright Laboratory Control Plus Freezer 232L ATEX

Gross Capacity (litres)	232L
Net Capacity (litres)	212L
Unit External Width (mm)	595mm
Unit External Depth (mm)	634mm
Unit External Height (mm)	1570mm
Product Weight (kg)	60kg
Packed Width (mm)	657mm
Packed Depth (mm)	734mm
Packed Height (mm)	1585mm
Packed Weight (kg)	65kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	480mm W x 1382mm H x 445mm D (step 175mm D x 255mm H)
Temperature Range	-18°C to -25°C
Colour	White/green
External Finish	White
Internal Finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, reversible
Door Open Depth (mm)	1185mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 front, 2 rear rollers
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	4 drawers, 2 fixed flaps
Number of Shelf Positions	Fixed
Shelf Size (mm) (WxDxH)	N/A
Shelf Material	Plastic
Number of Baskets	0
Basket Material	N/A
Basket Size (mm) (WxDxH)	N/A
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	Manual
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

## Specification: Upright Laboratory Control Plus Freezer 232L ATEX

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	<45dB
Energy Consumption (kwh)	-
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	60mm
Insulation Type	Cyclopentane
Blowing Agent	Yes
Anti Microbial Handle	No
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.