



Upright Laboratory Control Plus Fridge-Freezer Combi 263L ATEX

GDPA code: LSFC263UK-ATEX Lab ATEX UP F/F

SKU code: 444410677

- ✓ Operating Temperature Range of +2°C to +10°C for Fridge and -12°C to -20°C for Freezer
- ✓ Spark Free
- ✓ Advanced Controller
- ✓ Unique Anti Microbial Handle
- ✓ Audible Warning Alarms



The Experts in
Medical Refrigeration

Specification: Upright Laboratory Control Plus Fridge-Freezer Combi 263L ATEX

Gross Capacity (litres)	263 (Fridge: 159, Freezer: 104)
Net Capacity (litres)	245L (fridge 157L, freezer 88L)
Unit External Width (mm)	540mm
Unit External Depth (mm)	600mm
Unit External Height (mm)	1680mm
Product Weight (kg)	58kg
Packed Width (mm)	597mm
Packed Depth (mm)	684mm
Packed Height (mm)	1700mm
Packed Weight (kg)	64kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	Fridge 457mm W x 725mm H x 450mm D (step 220mm H x 170mm D), Freezer
Temperature Range	Fridge 2°C to 10°C, freezer -18°C to -25°C
Colour	White/green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1110mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	2 fitted locks (one on fridge, one on freezer)
Adjustable Feet (number & position)	3 adjustable front feet, 2 fixed rear
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	3 (Fridge)
Number of Shelf Positions	6
Shelf Size (mm) (WxDxH)	308mm W x 4.5mm H x 448mm D
Shelf Material	Wire
Number of Baskets	2 Solid half (Fridge)
Basket Material	2 full + 2 half
Basket Size (mm) (WxDxH)	-
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	Fridge: Automatic, Freezer: Manual
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

Specification: Upright Laboratory Control Plus Fridge-Freezer Combi 263L ATEX

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	<45dB
Energy Consumption (kwh)	0.91 kWh/24hours
Voltage Supply	230
Aux Output	230
Insulation Thickness (mm)	230
Insulation Type	Cyclopentane
Blowing Agent	Yes
Anti Microbial Handle	No
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.