



British healthcare and laboratory  
refrigeration specialists



## Laboratory Range

Under Counter Fridge and Freezer



FEATURES	KEY BENEFITS
Spark free interior	Our range of laboratory refrigerators and freezers offers the safe and reliable storage of volatile liquids, reagents, chemicals, poisons and samples with a spark-free interior that reduces the risk of internal explosion.
Lock	Security is essential in laboratories so our Lec Medical fridge and freezer come complete with door lock and two keys, ensuring that only authorised individuals are able to access them.
12 month warranty	Our under counter fridge and freezer come complete with a 12 month warranty for added peace of mind.




British healthcare and laboratory  
refrigeration specialists

## Laboratory Range

Under Counter Fridge and Freezer



### SPECIFICATION TABLE



Product Description		Laboratory Fridge
Product Dimensions (external)	mm	545 x 554 x 845mm (W x D x H)
Capacity	litres	130
Colour/Finish		White
Reversible Door	Y/N	Yes
Door Lock	Y/N	Yes
Drawer Number		1 Plastic
Shelf Number		3 Glass
Handle Type		Recessed
Defrost Capability		Auto Defrost
Annual Energy Consumption	kwh	117
Thermostat		Adjustable
Spark Free	Y/N	Yes
HFC Free	Y/N	Yes
CFC Free	Y/N	Yes
Product Description		Laboratory Freezer
Product Dimensions (external)	mm	545 x 554 x 845mm (W x D x H)
Capacity	litres	85
Colour/Finish		White
Reversible Door	Y/N	Yes
Door Lock	Y/N	Yes
Drawer Number		3 Plastic
Shelf Number		1 Plastic
Handle Type		Recessed
Annual Energy Consumption	kwh	157
Thermostat		Adjustable
Spark Free	Y/N	Yes
HFC Free	Y/N	Yes
CFC Free	Y/N	Yes



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.