

GG 4010

Gross capacity, total	382 l
Product dimensions (H/W/D)	190 / 59.7 / 68.0 cm
Energy consumption in 24h	1.176 kWh / 24h
Insertion slot	—



General product information

Output and consumption

Energy consumption in 24h	1.176 kWh / 24h
Net capacity, total	348 l
Freezer compartment	348 l
Gross capacity, total	382 l
Freezer compartment	382 l
Refrigerant	R 290
Cooling system, freezer compartment	static
Defrosting, freezing compartment	manual
Voltage	220-240 V ~
Frequency	50 Hz

## Dimensions and weight

Insertion slot	—
Product dimensions (H/W/D)	190 / 59.7 / 68.0 cm
Weight (with packaging)	88.0 kg
Weight (without packaging)	81.0 kg

## Control and functions

Control	Electronic control system
Control unit	Buttons
Temperature range	-14°C to -28°C
Temperature display	external digital
Malfunction: Warning signal	optical and acoustical
potential-free contact	—

## Design and features

Side wall material	steel
Colour	white
Door/Cover material	steel
Lock	fitted
Door hinges	right interchangeable
Material of interior containers	Plastic white
Storage shelves	7
Shelf material	Evaporator plates
Load capacity storage shelves, freezer compartment	60 kg
Door design	SwingDesign
Handle	Ergonomic slimline handle
Material of adjustable feet	Stainless steel



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.