

LGPv 1420

MediLine



Total gross capacity	1427 litre
Ext. dimensions in mm (w/d/h)	1430/830/2150
Int. dimensions in mm (w/d/h)	1250/700/1550
Energy efficiency class	
Energieverbrauch in 365 Tagen	3244 kWh
Climate rating ¹	5
Cooling system	–
Defrosting	automatic
Temperature range	-10°C to -26°C
Gastronorm	- mm
Housing material/colour	
Door/lid material	Steel
Liner Material	304-grade stainless steel
Type of control	electronic with text display
Temperature display	external digital
Interior light	–
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	8
Adjustable Shelving	–
Number of drawers	–
Number of baskets	–
Maximum number of baskets	–
Shelf material	grid shelves plastic-coated
Shelf loading	60 kg
Height of adjustable legs	-
Castors/type of castors	yes/swivel castors
Door handle	full-length integrated door handle
Foot pedal opener	–
Self-closing door	yes
Insulation	83 - 60 mm
Defrost water drain	–
Refrigerant	R 290
Airborne acoustical noise emissions	60 dB(A)
Gross/net weight	262,0/219,0 kg
EAN-Nr.	9005382137215

Subject to modification. Timestamp: 03.11.2016

LIEBHERR



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.