

# LKexv 3910

## MediLine



<b>Total gross capacity</b>	360 litre
<b>Ext. dimensions</b> in mm (w/d/h)	600/615/1840
<b>Int. dimensions</b> in mm (w/d/h)	440/435/1635
<b>Energy consumption</b> in 24 h	0,865 kWh
Climate rating <sup>1</sup>	SN-ST
Cooling system	—
Defrosting	automatic
Temperature range	+3°C to +16°C
Gastronorm	—
Housing material/colour	
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Interior light	—
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	6
Adjustable Shelving	5
Number of drawers	—
Number of baskets	—
Maximum number of baskets	—
Shelf material	glass
Shelf loading	40 kg
Height of adjustable legs	-
Castors/type of castors	—
Door handle	Ergonomic door handle
Foot pedal opener	—
Self-closing door	yes
Insulation	40 - 40 mm
Defrost water drain	—
Refrigerant	R 600a
Airborne acoustical noise emissions	48 dB(A)
Gross/net weight	74,0/68,0 kg
EAN-Nr.	9005382181898

Subject to modification. Timestamp: 26.04.2016

# LIEBHERR



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.