

# Ductaire Ducted Fume Cupboard



## Features & Benefits

**INFO**

- Single Walled Construction
- Air Bypass / Constant Air Volume
- BS : EN 14175 Compliant
- Auto Airflow Monitoring
- Digital Control Interface
- Custom Build Specifications

## High Specification Fume Extraction

### The Solution

The Ductaire Ducted Fume Cupboard offers robust, single-wall steel construction and advanced airflow control to deliver safe and efficient fume extraction.

Designed with a spacious internal working area and intuitive user interface, it ensures operator protection and ease of use.

Its streamlined layout and energy-efficient features help meet stringent safety standards while maintaining cost-effective operation.

### How it Operates

An Air Bypass / Constant Air Volume [CAV] system guarantees full operator protection against hazardous substances.

Negative pressure, sustained by a remote exhaust fan, ensures a constant, safe inflow of air inside the working chamber.



## Technical Information

**INFO**

	D1000	D1200	D1500	D1800
<b>External Dimensions*</b> (W x D x H)	1000 x 760 x 1300mm	1200 x 760 x 1300mm	1500 x 760 x 1300mm	1800 x 760 x 1300mm
<b>Internal Dimensions</b> (W x D x H)	950 x 550 x 1050mm	1150 x 550 x 1050mm	1450 x 550 x 1050mm	1750 x 550 x 1050mm
<b>Air Volume</b>	610m <sup>3</sup> /hr	750m <sup>3</sup> /hr	965m <sup>3</sup> /hr	1180m <sup>3</sup> /hr
<b>Lighting</b>	2no. 8W 2ft LED	2no. 8W 2ft LED	2no. 8W 2ft LED	2no. 8W 2ft LED
<b>Duct Spigot</b>	250mm Ø	250mm Ø	250mm Ø	250mm Ø



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.