

Chest freezer LT C200



- Low cost of ownership
- Limited maintenance
- On-site service

- Easily installable & movable
- Efficient cyclopentane insulation
- ISO 9001, 14001 & 45001 Certified

| Specifications | |
|--|-----------------------------|
| Technical | |
| Gross volume | 198 litres |
| Net volume | 189 litres |
| Temp. range | -25°C to -45°C |
| Temp. probes | 1/PTC |
| Compressor | Embraco |
| Baskets (incl./optional) | 1/6 |
| Refrigerant | R 290 |
| Noise level | 55 dB |
| Power consumption 24/h | 2,8 kWh |
| Power consumption | 1.022 kWh / year |
| Voltage | 230V |
| Frequency | 50 Hz |
| Features | |
| Audible & visual alarm | YES |
| Castors | 4 |
| Integrated lock | YES |
| Inner lids | NO |
| Defrost | Manual |
| Battery Backup | NO |
| USB | NO |
| Data connection | NO |
| Digital controller with display | XR30CX |
| Data logging | NO |
| High/Low temperature alarm | YES |
| Porthole for external probe & CO2 backup | Ø12 mm |
| Connection for remote alarm | Optional |
| Material | |
| Inner | Precoated white steel |
| Outer | Galvanized pre-coated steel |
| Dimensions | |
| Height (with castors) | 885 mm |
| Width | 923 mm |
| Depth (incl. condenser) | 630 mm |
| Depth (incl. condenser & handle) | 680 mm |
| Net weight | 63 kg |
| Internal dimensions | |
| Height | 635/385 mm |
| Width | 755/550 mm |
| Depth | 440 mm |
| Insulation | 80 mm |
| Insulation type | PUR |
| Internal humidity | <5 % |

LT C200

Product no. N111006

| Storage capacity | | | | |
|------------------|--------------|--------------|--------------------|-------------|
| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
| 50 mm | 6 | 3 | 18 | N311200 |
| | 11 | 9 | 99 | N311220 |
| | | | = 117 | |
| 75 mm | 4 | 3 | 12 | N311300 |
| | 8 | 9 | 72 | N311320 |
| | | | = 84 | |
| 100 mm | 3 | 3 | 9 | N311400 |
| | 5 | 9 | 45 | N311420 |
| | | | = 54 | |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.