



Model:

CL-01

Typology:

ILLUMINATOR

Description:

Cold light generator based on 100W dichroic halogen lamp.

Light source:

Dichroic halogen lamp 12V 100W, socket 6.35 (model OSRAM HLX 64627 or equivalent).

Controls:

ON/OFF switch and brightness adjustment on external controller.

Power Supply:

Input: 230Vac 50/60 Hz (optional: input 110Vac). Fuse 2AT

Accessories:

The output can be coupled to the following accessories:

- **CL-11:** double-arm optical fiber guide

Arm length 500 mm. Fitted with a fixed condensing lens at the end of each arm.

- **CL-12:** ring optical fiber guide

Guide length 700 mm. Guide diameter 16 mm. Diameter of the fixing ring: 55mm.

The circular end is suitable for all stereo series.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.