



Model:

**CL-31**

Typology:

**ILLUMINATOR**

**Description:**

*Double flexible arm (L = 60 cm) with double LED 2W.*

**Chromaticity coordinates (x,y):**

0.320, 0.320 (pure white)

**Optical power:**

Luminous flux: 160 lumen

Illuminance (at 10 cm distance): 50000 lux

**Aperture angle:**

35°

**Power supply:**

Brightness control potentiometer.

External power supply provided.

Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.