



Model:
CL-16.1

Typology:
ILLUMINATOR

Description:
LED ring light illuminator.

Cromaticity coordinates (x,y):

0.320, 0.320 (pure white), 6300K

Optical power:

Luminous flux: 500 lumen
Illuminance (at 10 cm distance): 6.000 lux (on 10cm diameter).

Power supply:

Knob for continuous brightness adjustment.
Different illumination settings: omnidirectional (ring completely on), oblique (only a quarter or half ring switched on).

External power supply provided, with built-in illumination controls.
Input: 100/240 Vac 50/60 Hz; output: 12 Vdc 150 mA

Accessories:

Fixing ring for stereo heads.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.