

*144-LED ring illuminator, EU plug*

Chromaticity and LED	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	LED number	144
	Brightness control	Manual
	Quadrant ON/OFF control	Yes
	Lifetime (hours)	> 65,000
	Max. required power (W)	2
Optical power	Luminous flux (lumen)	120
	Illuminance (at 10 cm distance, lux)	6,000 (on 10 cm diameter)
Power Supply for Illumination	Type	Internal
	Power plug type	EU plug
	Input voltage	100/240 Vac, 50/60 Hz, 150mA
	Output voltage	10-12Vdc
Accessories Included	Fixing ring	Yes



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.