

# LAB

Seat or stool options for use in laboratory settings that are easily cleaned for maintaining hygiene



#### BACK ANGLE ADJUSTMENT

Lift the back right lever to free the tilt angle of chair back to your preferred position.  
To lock in position push the lever back down

#### SEAT HEIGHT ADJUSTMENT

Lift the front right lever to operate the gas lift.  
This will increase the chair height.  
Lifting whilst seated will lower chair height

## DIMENSIONS

PLAB(1)  
Black polyurethane lab chair

SEAT HEIGHT 430-570mm  
SEAT WIDTH 460mm  
SEAT DEPTH 430mm  
OVERALL HEIGHT 830-960mm  
OVERALL WIDTH 630mm  
OVERALL DEPTH 630mm

PLAB(2)  
Black polyurethane lab chair

SEAT HEIGHT 430-570mm  
SEAT WIDTH 460mm  
SEAT DEPTH 430mm  
OVERALL HEIGHT 830-960mm  
OVERALL WIDTH 630mm  
OVERALL DEPTH 630mm

PLAB(1)D  
Black polyurethane lab chair  
with draughtsman base

SEAT HEIGHT 540-800mm  
SEAT WIDTH 460mm  
SEAT DEPTH 430mm  
OVERALL HEIGHT 940-1200mm  
OVERALL WIDTH 630mm  
OVERALL DEPTH 630mm

PLAB STOOL  
Black polyurethane lab stool

SEAT HEIGHT 430-570mm  
SEAT WIDTH 450mm  
SEAT DEPTH 410mm  
OVERALL HEIGHT 530-670mm  
OVERALL WIDTH 630mm  
OVERALL DEPTH 630mm

PLAB SMALL STOOL  
Black polyurethane small lab stool

SEAT HEIGHT 410-540mm  
SEAT WIDTH 350mm  
SEAT DEPTH 350mm  
OVERALL HEIGHT 410-540mm  
OVERALL WIDTH 630mm  
OVERALL DEPTH 630mm



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.