


**1. Part Number:** A836-250

**2. Photograph:**



### 3. Product Information:

Description:	Stainless Bottle 250 ml
Material of Construction:	Body: 316L stainless steel Lid: 316L stainless steel Clamp: 304 stainless steel Gasket: Silicone (conforms to FDA CFR 177.2600)
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	250ml
Brimful Volume:	380ml
Dimensions:	A = 64 mm B = 48 mm C = 97 mm D = 83 mm
Nominal Weight (body, lid, clamp & gasket):	710 g
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

Note: The neck of the container is fitted with a 2" diameter ferrule which conforms to BS4825-3

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Feb-21



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.