

**1. Part Number:** A200-200

**2. Photograph:**



**3. Product Information:**

Description:	Pallet Knife with Plastic Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Black PVC
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	200mm
Blade Width:	30mm
Overall Length:	310mm
Nominal Weight:	59g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Apr-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.