



**1. Part Number:** A314-100

**2. Photograph:**



### 3. Product Information:

Description:	Ice Scoop
Material of Construction:	304 Stainless Steel
Method of Construction:	Spot Welded
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	180mm
Body Length:	100mm
Body Diameter:	50mm
Diameter of Drainage Holes in scoop body:	3mm
Nominal Weight:	90g
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.