

**1. Part Number:** A376-500

**2. Photograph:**



**3. Product Information:**

Description	Weighing Scoop 500ml
Nominal Volume:	500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	320mm
Body Length:	205mm
Maximum Width of Scoop Body:	110mm
Minimum Width of Scoop Body:	65mm
Nominal Weight:	347g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.