



1. Part Number: A534-1000

2. Photograph:



3. Product Information:

Description	Stainless Mixing Bowl 1000ml
Nominal Volume:	1000ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece construction
Surface Finish:	Better than 0.5 microns Ra
Height	86mm
Top Diameter:	170mm
Weight:	140g
Graduation marks:	None
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.