

**1. Part Numbers:** A586-12

**2. Photograph:**




### 3. Product Information:

Description	Heavy Duty Bucket 12L
Nominal Volume:	12L
Brimful Volume:	12L
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior Lugs for handle are spot welded to exterior of bucket
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	310mm
Bucket Height:	262mm
Graduations:	No
Skirt Band:	No
Nominal Weight:	2.4 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
 Position: Quality Manager  
 Date: 26-Oct-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.