



	Product Specification	Specification Number	A686-15L
		Revision	01

Part No. A686-15L
Description Lid for 15L Bucket with Spout
Material of Construction: 316L stainless steel



Method of Construction: Handle spot welded to top of lid
Crevice free interior
Surface Finish: Better than 1.0 microns Ra
Max. width over spout: 380mm
Diameter: 330mm
Nominal Weight: 1kg
Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 6-Aug-19
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.