

1. Part Number: A716-300

2. Photograph:



3. Product Information:

Description:	Powder Funnel 300mm
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	8 Litres
Top Diameter:	300mm
Top Diameter including flange:	330mm
Overall Height:	310mm
Spout Length:	25mm
Spout Diameter:	Spout ends with a 2" ferrule (BS4825-3)
Nominal Weight:	1720g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Nov-19



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.